

# MaggotTracker

*Drosophila melanogaster* larval locomotion

## User manual for MaggotTracker software components.

### MagAnalyzer.

#### I. Overview

MagAnalyzer calculates locomotion parameters of one or more folders containing videos previously processed by MagRecognizer.

#### II. Installation

Download MagAnalyzer ([downloads](#)).

Simply unzip the ZIP file. The Java programs are compiled already.

#### III. Input files

The input is one folder name, normally containing many sub-folders that were already processed by MagRecognizer.

Any sub-folder with name that starts with 'bad' is skipped.

#### IV. Running the software

There is no user-interaction when MagAnalyzer processes a folder containing videos. Hence, the method we use to run it is via command-line:

On *Windows* systems:

- open a "Command Prompt" window
- Go to the folder where MagAnalyzer is located. For example: `cd \maggottracker\maganalyzer`
- Type 'MagAnalyzer' followed by the folder containing the (processed) videos. For example: `run c:\data\processed`

On *Mac OSX*, or *Linux* systems:

- open a "Terminal" window
- Go to the folder where MagAnalyzer is located. For example: `cd /Users/ana/maggottracker/maganalyzer`
- Type './linux\_MagAnalyzer.sh' followed by the folder containing the (processed) videos. For example: `./osx_run.sh /Users/ana/data/processed`

#### V. Result file

A single plain-text tab-separated file named datadm.txt is created with a header line listing all parameter names, and one line of parameter values for each video.

The file format facilitates loading data into a database, or into a spreadsheet.

---

[Home](#) | [Download](#) | [Documentation](#)